




Smart-Mini Water Purification System Operation Manual



Version:1612

Notice

If it is the first time for you to use our product, please carefully read this manual which will give you a lot of help. We take responsibility for regular maintenance and repair work instead of consequences caused by improper operation.

The content of the publication is subject to change and/or updating without notice. The content of this manual is only for reference use. The manufactory and the distributors shall not be responsible for losses caused by incorrect description or misuse of this manual. The graphics contained in the manual are as general representative which may have some errors compared to entity. But the description is definitely match the function.

In accordance with copyright law, the company holds and retains the exclusive right to all the works. Without written consent provided by company, any other organization or individual has no right to any form of manual changes.

Preface

Thank you for using the water purification system designed and manufactured by Shanghai Canrex Analytic Instrument Co, Ltd.

If you have any good suggestion, Please contact us, we will improve our products and after sale service continuously.

Global Exclusive Agent: Healforce development (hongkong) Co. ltd. Tel: (00852) 28987303

China Exclusive Agent: Nison Instrument (Shanghai) Co. ltd. Tel: (8621)62728646

Manufacture: Shanghai Canrex Analytic Instrument Co. ltd. Tel: (8621)50911997

ISO 9001:2008 Certificate

TUV NORD

CERTIFICATE

Management system as per GB/T 19001-2008/ ISO 9001 : 2008

In accordance with TUV NORD CERT procedures, it is hereby certified that

SHANGHAI CANREX ANALYTIC INSTRUMENT CO., LTD.
3 Floor, 4 Factory Building, No. 298, Lianzhen Road, Pudong New
District, Shanghai, P. R. China

Organization Code: 78112153-3

applies a management system in line with the above standard for the following scope

Design and Manufacturing of Water Purification Systems for Laboratories

Certificate Registration No. 04 100 062186

Audit Report No. 2.5-4951/2015

Valid from 2015-06-12

Valid until 2018-06-11

Initial Certification 2006-06-23



Certification Body
at TUV NORD CERT GmbH

Printing Address
TUV NORD Hangzhou Co., Ltd
Shanghai Branch
11F, Sail Tower, No 266 Han Kou Rd.,
Huang Pu District, Shanghai, China
200001
Issue Date: 2015-04-29

This certification was conducted in accordance with the TUV NORD CERT auditing and certification procedures and is subject to regular surveillance audits.

TUV NORD CERT GmbH

Langemarckstrasse 20

45141 Essen

www.tuev-nord-cert.com



Note: This certificate information can be queried at the official Website of Certification and Accreditation Administration of the People's Republic of China: www.cnca.gov.cn

SH0_F04001_PG_0182_Sector 2

Contents

Chapter 1 General	1
1.1 Safety information	1
1.2 System operating environment	1
Chapter 2 System Introduction	2
2.1 Specifications for the system	2
2.2 Diagram	2
Chapter 3 Instructions before installation.....	3
3.1 check before installation.....	3
3.2 Option spare parts.....	3
3.2.1 Pressure limiter CR-SP829.....	3
3.2.2 Enhanced pretreatment cartridge CR-SP102L	3
3.2.3 Remote water dispenser CR-SPRWD-1	3
Chapter 4 Installation	4
4.1 Connection of feed water to the system	4
4.2 Connection of the drain tubing	4
4.3 Connection of the overflow tubing.....	4
4.4 Installation of the pretreatment cartridge.....	4
4.5 Installation of the purification cartridge	5
4.6 Installation of the ultra-purification cartridge	5
4.7 Installation of the final filter	6
4.8 Installation of the leakage protection sensor	7
4.9 System electrical connection	7
Chapter 5 Parameter setting and operation.....	8
5.1 Control panel area.....	8
5.2 LCD display area	8
5.3 symbols and graphics for LCD display	8
5.4 Function Keypad	9
5.5 Operation	9
5.6 Parameters setting.....	10
5.6.1 Alarm Settings	11
5.6.2 Volume of water setting	13
5.6.3 Calendar setting	14
5.6.4 Ultra-pure water quality setting.....	14
5.6.5 Unit selection.....	15
5.6.6 Dispensing flow setting	15
Chapter 6 Troubleshooting	17
6.1 No display.....	17
6.2 Making RO water is too low	17
6.3 The pure water conductivity is too high	17
6.4 Ultra-pure water resistivity is too low	17
6.5 The flow rate of ultra-pure water is too low	18
Chapter 7 Order information of option parts.....	18
Chapter 8 Appendix	19

Appendix 1 Fuse replacement	19
Appendix 2 pressure reductor installation	20
Appendix 3 inlet valve unit installation.....	21
Appendix 4 Installation of the remote water dispenser	22
Appendix 5 Components connection to main-board	23

Chapter 1 General

1.1 Safety information

You must use the safety norms according to this manual before using the Smart-Mini system, especially water and power supply. It is necessary to refer to this manual when you install or operate the Smart-Mini system. Unqualified using environment will endanger the normal operation, or even damage the whole system.

The installation, commissioning and maintenance of the equipment can only be completed by Heal Force or its authorized agents. Heal Force did not responsible for duties and responsibilities if the equipment is disassembled by unauthorized dealer or service personnel.

1.2 System operating environment

- Indoor sue
- Avoid direct sunlight
- Between 0-2000 meters altitude
- Ambient temperature: 4-45°C
- Operating Voltage: 220-230V AC 50/60Hz
- Inlet pressure: 0.1-0.4MPA.
- The fluctuation range of the main power supply can be 10% of the normal voltage.
- Transient over-voltage power supply is grade 2.
- There is an drains around the equipment(1 meters)
- Adequate indoor ventilation
- Mounting surface must be fixed, level, and not burning
- Avoid direct sunlight
- No heat source next to the equipment.
- Be away from strong magnetic field.

Chapter 2 System Introduction

2.1 Specifications for the system

Feed water requirement	Conductivity $\mu\text{S}/\text{cm}@25^\circ\text{C}$	< 400
	Pressure MPa	0.1~0.4
	Temperature $^\circ\text{C}$	5~40
Pure water	Ion rejection rate $\%@25^\circ\text{C}$	≥ 95
	Bacteria remove %	> 99
	Make rate* L/H@25 $^\circ\text{C}$	15
Ultra-Pure water	resistivity $\text{M}\Omega.\text{cm}$	> 18
	TOC** ppb	1-10
	microorganisms*** (CFU/ml)	≤ 1
	Particles**** (0.22 μm) /ml	< 1
	Dispensing rate L/H	60-90
	Power Consumption W	120
	Weight kg	25
	Storage tank	6L
	Main System Dime. W*D*H	525*314*571

Note:


* The making rate will be changed with temperature variation. (3% per 1 $^\circ\text{C}$)

** The feed water TOC is less than 1000ppb.

*** with 0.22 μm final filter cartridge

**** with Ultra-filter cartridge

2.2 Diagram



- | | | | | | |
|----|------------------------|----|------------------------------|----|--------------------|
| 1 | pre-filter | 2 | low pressure sensor | 3 | inlet valve |
| 4 | conductivity sensor 1 | 5 | pretreatment cartridge | 6 | RO pump |
| 7 | RO membrane | 8 | conductivity sensor 2 | 9 | water tank |
| 10 | dispensing pump | 11 | high pressure sensor | 12 | UV lamp |
| 13 | purification cartridge | 14 | ultra-purification cartridge | 15 | resistivity sensor |
| 16 | micro-filter | 17 | dispensing valve | 18 | dispenser(option) |
| 19 | RO flushing valve | 20 | wastewater valve | 21 | drain valve |

Chapter 3 Instructions before installation

3.1 check before installation

Users please ensure whether the product you received is your purchase before installation.

You had better check all items according the packing list.

- Whether the power site is close to the system.
- Whether the feed water meets the requirements.
- Whether the feed water pressure is normal.
- Whether the feed water supply is nonstop.
- Whether the feed water pipe has 1/2" NPTF whorl connector.
- Whether the strength of the wall can bear the system (including tank full of water) weight.
- Whether the system's height can clearly observe both LCD screen and keypad operations on control panel.
- Whether there is adequate space to substitute supplies and connecting pipe as well around and behind the system.
- Whether you have installed more than 1/2" drains, not more than one meter far away from the system, intended for discharging waste-water and water overflowing.
- Whether the installment circumstances are consistent with the regulations.
- Whether there is reserved space if you are looking forward to having remote water dispenser and pedal switch installed.

3.2 Option spare parts

3.2.1 Pressure limiter CR-SP829

Please purchase this product additionally when the feed water pressure is too high ($> 0.4\text{MPa}$) . This component can reduce pressure to 0.4Mpa below.

3.2.2 Enhanced pretreatment cartridge CR-SP102L

Enhanced pretreatment cartridge is a multi-stage pretreatment module to removes the impurities in tap water including suspended particulates, colloids, microorganisms, organics, chlorine and heavy metals to ensure the following purification modules work properly.

3.2.3 Remote water dispenser CR-SPRWD-1


For Multiple water dispensing, the longest distance is 3 meters.

Chapter 4 Installation

It is important for the user to master installation and removal of quick connector. If user does not operate properly will lead to water leakage, quick connector damage and pipeline broken.

4.1 Connection of feed water to the system

Before connecting, close any valves supplying feedwater to the system.



A 1/2" male fitting

D pre-filter


B 1/2" to 1/4" connector

E 1/4" blue tubing(2m)

C 1/4" blue tubing(20cm)

F "FEED"

4.2 Connection of the drain tubing



A "DRAIN1"


D "DRAIN2"

B 1/4" black tubing(2m)

E 1/4" black tubing(2m)

C floor drain

4.3 Connection of the overflow tubing




A "OUT"

B 1/4" black tubing(2m)

C floor drain

Note: The overflow tubing must have a certain slope. It shall not be bent or arched.


4.4 Installation of the pretreatment cartridge



Step A: unscrew the screw of the side door with a screwdriver.

Step B: pull the handle of the side door with both hands, forced back to pull the a side door out of 10mm ~ 13mm.

Step C: carry out the side door.



Pretreatment cartridge is installed in the right chamber of the system.

Get the pretreatment cartridge from the box (make sure number is correct for the CR-SP101M), open the vacuum bag, remove protection cover for the inlet and outlet of pretreatment cartridge. Make sure there is a black "O" ring in the inlet and outlet, which is in the bottom of the hole.

Place the square tenon at the low of the pretreatment into the corresponding square hole on the main unit, shown as in the figure "D" above.

Push the joint at the up of the pretreatment cartridge to the main unit shown as in the figure "E" above. Make sure the "O"-rings of pretreatment cartridge can be completely sealed with the joints.

There is locking plate connected with nylon line on the joints of main unit. Insert the locking plate into the slot of metal bar, as shown in the figure F above.

4.5 Installation of the purification cartridge


Remove the left side door (front view).

Purification cartridge is installed in the left chamber of the system.


Get the purification cartridge from the box (make sure number is correct for the CR-SP301M).

The operation method: please refer to section 4.4 of this chapter.

4.6 Installation of the ultra-purification cartridge



Press and pull the latch at the top of the front door, to remove the front door (press down about 3mm , do not press down hardly which may cause deformation or even fracture bayonet). As shown in figure "A""B"



Left the upper and lower latches, As shown in Figure C.


Take out the ultra-purification cartridge from the box (make sure number is correct for the CR-SP303M), open the vacuum bag, remove protection cover for the port of cartridge. Make sure there is a black "O" ring in the ports, which is in the right position.(should be close to the inlet and outlet inside the convex shoulder). If the position is not right, Please adjust it to the right position.

Push the ultra-purification cartridge firmly into place , Both the top and bottom ports of the cartridge have to be fully seated shown as figure "D" above.

Place the latch as origin. Make sure that both the top latch and bottom latch are fully closed. Each latch should "click" when it is fully closed.

Install the front door as this chapter. The purification cartridge installation is completed

4.7 Installation of the final filter




Take the final filter and quick connector (with 40mm 1 / 4 "white hose) from the attachment bag.

Wind sealing tape around external screw for 2-3 laps then screw on quick connect of filters, and tighten. Shown in Fig A

Insert the white hose into the water producing terminal of the system.

4.8 Installation of the leakage protection sensor




Get the cable one end connect to hole at the back of main unit. Put the leakage protection sensor on the table near the main unit.

4.9 System electrical connection

Using the supplied electrical power cord, connect it to the system and then to an earth grounded outlet.

Chapter 5 Parameter setting and operation


5.1 Control panel area



A: dashed box for the LCD display area.

B: dashed box for the function keypad board.

5.2 LCD display area



A: system model and the current temperature display .

B: current status

C: product water quality and Volume dispense status.

D: date and feed water conductivity.



E: graphic for display status.


5.3 symbols and graphics for LCD display


FEED Feed water conductivity
 UV UV light
 Vol.Disp. Fixed Volume Dispensing
 PF Pretreatment cartridge


RO Pure water conductivity
 DS RO deionization rate
 UP Products resistivity
 MF Micro-filter


 Without feed water


 Feed water low quality

FE
 RO failure

RO
 UV light failure


UV
 RO failure


DS
 Ultrapure water recycle


 System leakage


 Pretreatment failure

PF
 Micro-filter failure

MF
 Purification cartridge failure

DI
 Ultrapure water producing

 Water tank level

5.4 Function Keypad

There are 4 function keypads in the system: DISP./▲, Flush/▼, Stop/↵ and Menu

- DISP./▲ Dispense product water keypad at the operation, it will be become adjusting keypad at the menu mode.
- Flush/▼ RO cartridge flush keypad at the operation, it will be become adjusting keypad at the menu mode.
- Stop/↵ Stopping water dispensing or flushing at the operation, it will be become confirmation keypad at the menu mode.
- Menu Parameter setting keypad, it will be become return keypad at the menu mode.


5.5 Operation

Turn on feed water supply. Ensure that all the pipe connecting parts have no leaking. Switch on the power supply at the back of the main unit.


After the main unit is powered on, LCD will display a welcome message, which is followed by the current state of the system. The system will remain in flushing status about 1 minute. As shown below:

Smart-Mini	25.0°C
RO member Flushing...	
1:00	
08:00 20/09/2016	

When the flushing process of main unit is finished, the system will automatically enter the water producing state. The system will flushing automatically RO membrane before and after producing RO water. The system also shows water conductivity and RO deionization rate. As shown below:

Smart-Mini	25.0°C
RO water making	
RO: 10.0 μ s/cm	
DS: 98%	
20/09/2016 FEED: 400 μ s/cm	


When the storage tank of the system is full , the system will enter the standby mode automatically. As shown below:

Smart-Mini	25.0°C
Standby	
Ideal for lab	
08:00 20/09/2016	


Press "DISP./▲" keypad you can take the pure water.



Ensure that the pure water at least 20% of the water tank.

Smart-Mini	25.0°C
Ultra-water Dispensing	
UP: 18.2 MΩ · cm Vol Disp: OFF	
08:00 20/09/2016	


If the quality of current ultrapure water can not meet the requirement of original setting (ie 16.0MΩ.cm), the system will automatically enter recycling mode to improve the water quality. And if the water quality reaches the set value, the main unit starts to produce ultra-pure water. As shown below:

Smart-Mini	25.0°C
Ultra-water Cycle	
UP: 15.2 MΩ · cm Vol Disp: OFF	
08:00 20/09/2016	

Press "Stop/↵" keypad to stopping take product water.
The system has manual flush function for RO membrane cartridge, press the "Flush /▼" to flush RO membrane cartridge.

Press the "Stop/↵" keypad to stop flushing operation.

When the system is in the standby state, the ultrapure water in the loop will automatically recycle for 5 minutes every 2 hours, which will guarantee the quality of ultrapure water to the setting requirement. As shown below:

Smart-Mini	25.0°C
Standby	
Cycling...	
08:00 20/09/2016	

5.6 Parameters setting

Press "Menu" keypad to enter the menu, language selection, as shown below:

语言	LANGUAGE
中文	
ENGLISH	

Press "DISP./▲" or "Flush/▼" keypad to select "ENGLISH". Press the "Stop/↵" keypad to confirm the setting and to the next level. As shown below:

Password: 85

Press or hold down the "DISP./▲" or "Flush/▼" keypad to enter the correct password. The default password is 85. Enter the password, press the "Stop/↵" keypad to confirm the entry of next menu level, as shown below. If the password input is error, the system will be directly reset.

Parameter setting
Service setting

Press the "DISP./▲" "Flush /▼" keypad to select "Parameter setting". Press the "Stop/↵" keypad to enter the parameter setting level

After entering the parameter setting menu, you can see the contents of the setting are: Alarm settings, Calendar, Prod. unit, Volume of water, Prod. Resis., Dispensing flow. As shown below:

Alarm setting		
Alarm setting	Calendar	Prod. unit
Volume of water	Prod. Resis.	Dispensing flow

5.6.1 Alarm Settings

After entering the parameter setting menu, press the "DISP./▲" or "Flush/▼" keypad to select "Alarm setting". Press the "Stop/↵" keypad to confirm the settings menu to enter the alarm category including various settings, as shown below:

PF life :	400h	500h 600h
PF life	MF life	RO
UV life	Feed	DS

Alarm time setting for pretreatment failure

Press "DISP./▲" or "Flush/▼" keypad to select "PF life". Press the "Stop/↵" keypad to confirms. Press "DISP./▲" or "Flush/▼" keypad to select one of the 400,500,600hour according to the tap water quality. Press the "Stop/↵" keypad to confirms the setting and switches to the next menu level.

Select alarm time UV lamp failure

LCD display is as follows:

UV life : 4000h 5000h 6000h		
PF life	MF life	RO
UV life	Feed	DS

In this menu level the maintenance interval for the UV oxidator of main system must be reset after exchanging the UV lamp. the manufacture setting is 4000hour for UV lamp life time, the user can press “Disp./▲, Flush/▼” keypad to select one of the 4000,5000,6000,hour according to the user desire. Press the “Stop/↵” keypad to confirms the setting and switches to the next menu level.

Micro-filter failure time alarm options

LCD display is as follows:

MF life : 4000h 5000h 6000h		
PF life	MF life	RO
UV life	Feed	DS

In this menu level the maintenance interval for the Micro-filter Cartridge life time of main system must be reset after exchanging the cartridge. the manufacture setting is 4000hour for Micro-filter Cartridge life time, the user can press “Disp./▲, Flush/▼” keypad to select one of the 4000,5000,6000 hour according to the user de sire. Press the “Stop/↵” keypad to confirms the setting and switches to the next menu level.

Feed water quality alarm setting

LCD display shown as follows:

Feed Water Alarm-H : 400 uS/cm		
PF life	MF life	RO
UV life	Feed	DS

The limit of the feed tap water conductivity can be adjusted by the user. If the adjusted limit will be overstepped it will decrease the life time of Pre-treatment cartridge, RO membrane and purification cartridge. the manufacture setting is 400us/cm, the user can press “Disp./▲, Flush/▼” keypad to adjusting according to the user desire. Press the “Stop/↵” keypad to confirms the setting and switches to the next menu level.



When setting the value, please refer to the actual water quality. Only to eliminate the alarm for feed water. If the feed water quality is poor, the overall performance and lifetime of each module will be affected.

Pure water quality alarm settings

LCD display shown as follows:

RO Alarm-H : 20.0 uS/cm		
PF life	MF life	RO
UV life	Feed	DS

The limit of the RO water conductivity can be adjusted by the user. if the adjusted limit will be overstepped it will decrease the life time of purification cartridge. the manufacture setting is 20us/cm, the user can press “Disp./▲, Flush/▼” keypad to adjusting according to the user desire. Press the “Stop/↵” keypad to confirms the setting and switches to the next menu level



If this value is set too large, then the performance and lifetime of the ultrapure water purification cartridge will be affected.

Alarm settings for reverse osmosis desalting rate

LCD display is as follows:

DS Alarm-L : 90 %		
PF life	MF life	RO
UV life	Feed	DS

The limit of the desalting rate of the RO membrane can be adjusted by the user. The display of the system shows the alarm message “Ds light” if the adjusted limit will be overstepped (manufacture’s setting is 90%). press “Disp./▲, Flush/▼”keypad to increases or decreases the value, press the Stop/↵” keypad to confirm the setting and switches to the next menu level.

This value should not be set too low, otherwise the pure water quality will be affected. If the quality of inlet water is too poor, Please set high value for higher pure water requirement. Value setting range: 50 to 99.

5.6.2 Volume of water setting

After entering the parameter setting menu, press “DISP./▲” or “Flush/▼” keypad to select “Volume of water”.

Volume of water		
Alarm setting	Calendar	Prod. unit
Volume of water	Prod. Resis.	Dispensing flow

Press the “Stop/↵” keypad to confirms.

Volume: OFF		
Alarm setting	Calendar	Prod. unit
Volume of water	Prod. Resis.	Dispensing flow

Press “DISP./▲” or “Flush/▼” keypad to set the volume. Press the “Stop/↵” keypad to confirms the setting and switches to the next menu level.

5.6.3 Calendar setting

After entering the parameter setting menu, press “Disp./▲” or “Flush/▼” keypad to select the "Calendar".

Calendar		
Alarm setting	Calendar	Prod. unit
Volume of water	Prod. Resis.	Dispensing flow

In this menu level the user want to change the timer in main system ,you can be press the “Disp./▲, Flush/▼” keypad to set the correct date and time according to the mention of the display. Press the “Stop/↵” keypad to confirms the setting and switches to the next menu level.

Year 2016		
Year 2016	Month 11	Date 20
Hour 12	Minute 22	

5.6.4 Ultra-pure water quality setting

After entering the parameter setting menu according to this chapter, press “Disp./▲” or “Flush/▼” keypad to select the "Prod. Resis." function.

Product resistivity		
Alarm setting	Calendar	Prod. unit
Volume of water	Prod. Resis.	Dispensing flow

Press “Stop/↵” keypad to enter the water quality parameters setting interface, as shown below:

Prod. Resis. : 16.0 MΩ · cm		
Alarm setting	Calendar	Prod. unit
Volume of water	Prod. Resis.	Dispensing flow

Press “Disp./▲, Flush/▼” keypad to adjusting according to the user desire. Press the “Stop/↵” keypad to confirms the setting and switches to the next menu level.

5.6.5 Unit selection

By setting this parameter, user can choose to the unit form of the ultrapure water, MΩ.cm or μS/cm.

After entering the parameter setting menu according to this chapter, press “Disp./▲” or “Flush/▼” keypad to select the "Prod. Unit." function.

Product unit		
Alarm setting	Calendar	Prod. unit
Volume of water	Prod. Resis.	Dispensing flow

Press "Stop / ↵" button to confirm the unit selection interface, as shown below:

Prod. Unit : MΩ · cm μ S/cm		
Alarm setting	Calendar	Prod. unit
Volume of water	Prod. Resis.	Dispensing flow

In this menu level the user can select the display unit of product water , the user can press “Disp./▲, Flush/▼” keypad to selecting one of “MΩ.cm” or “us/cm” according to the user desire. Press the “Stop/Enter” keypad to confirms the setting and switches to the next menu level.

5.6.6 Dispensing flow setting

By setting this parameter, user can calibration the accuracy of fixed quantity of taking water such as 1L or 5L.

After entering the parameter setting menu ,Press the “Disp./▲” or “Flush/▼” keypad to select the "Dispensing flow" function.

Dispensing flow		
Alarm setting	Calendar	Prod. unit
Volume of water	Prod. Resis.	Dispensing flow

Press “Stop/↵” keypad to enter the water dispensing flow parameters setting interface, as shown below:

Dispensing flow : 1.20 L/m		
Alarm setting	Calendar	Prod. unit
Volume of water	Prod. Resis.	Dispensing flow

Press the "DISP/▲" "Flush /▼" keypad to enter the actual value.
 Press the "Stop/┘" keypad to confirm the setting and switch to the next menu level.
 Note: Suggestion the user had better to calibrating fixed quantity every half month.

Chapter 6 Troubleshooting

6.1 No display

Description	Solution
No power supply for power outlet	Confirm power outlet is normal
Electrical outlet loose	Change electrical outlet
System power cord is not properly installed	Make sure the power cord is plugged into the electrical outlet of main unit
system power is turned off	Turn on the power switch on the back of the main unit
Power cord is damaged	Change the power cord
Fuse burned	Change fuse
Internal component failure	Contact Technical service

6.2 Making RO water is too low

Description	Solution
No feed water supply, warning for water shut off	confirm the patency of external water supply
Pre-filter blocked, warning for no water	Rinse the filter and replace it when necessary
Feed water pressure is too low, indicating no water warning	Needs another pressure pump
water supply connection error, indicating no water warning	Make sure that feed water pipe connecting to the FEED connector of the main unit.
Pretreatment blocked	Replace the pretreatment module
There are impurities on the RO surface	Press "Flush/▼"button to flush RO
RO blocked	Replace the RO module
Pure water outlet connection error	Verify that the RO OUT of main unit connected with RO IN of water tank
Water temperature is too low	Add the preheat system
Internal component failure	Contact Technical service

6.3 The pure water conductivity is too high

Description	solution
Feed water conductivity is high	Choose RO-2 module
RO exhausted, showing RO or DS alarm	Replace the RO membrane
New reverse osmosis membrane with remaining protection liquid.	Press "Flush/▼"button to rinse the membrane
Internal component failure	Contact Technical service

6.4 Ultra-pure water resistivity is too low

Description	Solution
Purification cartridge is exhausted	Replace the Purification cartridge
The unit display conductivity value	Switch the unit mode to resistivity
Purification cartridge installation error	Make sure the purification cartridge installation OK.

Disinfectant is not fully discharged after disinfection	Discharge amount of the ultra-pure water
Internal component failure	Contact Technical service

6.5 The flow rate of ultra-pure water is too low

Description	Solution
System is in recycling mode	Wait for the producing water reaching preset quality
Water tank level is less than 20%	Wait for pure water level increase in water tank
Ultra-pure water resistivity is too low	Adjust the ultra-pure water quality to the requirement
Ultra-pure water resistivity value is setting too high	Adjust the ultra-pure water quality to the requirement
Supply pipe connection error	Confirm the RO OUT of water tank connecting to the RO IN of main unit
Internal component failure	Contact Technical service

Chapter 7 Order information of option parts

Description	specification	Order No.	Remark
Enhanced pretreatment device	Pretreatment cartridges	CR-SP102	
Pressure reductor	5KG	CR-SP829	
Remote water dispenser	Smart	CR -SPRWD-1	

Chapter 8 Appendix


Appendix 1 Fuse replacement

The fuse used by the system is 1A/250V, which is installed in the power outlet at the back of the main unit.

Pull the power cord from the external power outlet.

Unplug the power cord from the main unit. You can see a fuse holder below the power socket.

Use a screwdriver to pull the fuse holder from socket. As shown below




Take out the used fuse from the fuse holder, install new one into holder and push to the end.

Push the fuse holder to the square hole which lies in the lower part of power socket.

Fuse replacement is complete, plug the power cord and connect to an external power supply.

Appendix 2 pressure reductor installation

The feed water pressure for the system is normally 0.4MPA. If the pressure exceeds this value, a pressure reductor needs to be installed after the pre-filter, as shown below.



Please install the pressure reductor according to the water flow direction, do not install it in the wrong direction, otherwise it will not work normally.


The pressure reductor is adjustable, you can use tools to rotate the cross slots to change the feed water pressure relief. Please connect a pressure gauge to the outlet and block the following pipe when you adjust the pressure. Then you need to regulate slowly in accordance with the pressure gauge.

Do not rotate the regulator when you do not have any water pressure test instrument.

Appendix 3 inlet valve unit installation

Some of the laboratories have no threaded tap. The fast connection can not be used to connect to an external water supply pipe. The inlet valve components can be optional alternative, which is installed in the middle of external pipe.

Inlet valve subassembly is usually installed on the outside of the tap, as shown below.



Inlet valve subassembly includes a T-way and a ball valve. Installation steps are as follows

Turn off the external tap water supply valve.


Remove the tap from water supply pipe.

Take T-way which with teflon seal one end links to external water supply pipes, the second end to the tap and the third end to the ball valve with teflon seal.

Take the pre-filter connect to the ball valve.

Install other components in accordance with the manual

Appendix 4 Installation of the remote water dispenser



A Get a white tube (1 / 4 " , 3m), insert one end of the tube into the "RWD1-1" on the back of main unit, and the other end into the one of the connector on the back of RWD.


B Get a white tube (1 / 4 " , 3m), insert one end of the tube into the "RWD1-2" on the back of main unit, and the other end into the one of the connector on the back of RWD.

C Get a PS/2 cable, insert one end into the "RWD1" on the back of the system, and insert the other end into the PS/2 socket of the remote water dispenser.

Appendix 5 Components connection to main-board



The mainboard have 220V or 110V voltage. User should not arbitrarily open the main unit shelter.



Global Exclusive Agent

Nison Instrument (Shanghai) Limited

Address:

Floor 16th, 1065 West zhongsan Road, Shanghai, P. R. China

TEL: 86-21-62728646, 62712931

FAX: 86-21-62710529

E-mail: export@hf.healoo.com

Website: www.healforce.com

Manufacturer

Shanghai Canrex Analytic Instrument Co, Ltd.

Address:

No.4 Workshop Building, 298 Lian-Zhen Rd. Pudong Shanghai

P. R. China

