

INTRODUCING: THE MSOT SCANNER FAMILY.

Multispectral Optoacoustic Tomography (MSOT): The next generation in molecular imaging.



MSOT inSight64

MSOT inVision128

MSOT inVision256-TF

- Optoacoustic scanner for small animal imaging
- Anatomical, functional and molecular contrast
- Real-time image acquisition, 150 μ m resolution

MSOT: MULTISPECTRAL OPTOACOUSTIC TOMOGRAPHY

MSOT technology

Imaging sequence:

- Illumination of tissue with laser pulses at multiple wavelengths
- Detection of induced ultrasound pressure waves
- Spectral unmixing to analyze individual absorbers

Technology benefits:


Combines the molecular specificity of optical imaging with the depth and spatiotemporal resolution of ultrasound

• Molecular specificity:

Identify and quantify disease-related biomarkers, revealing endogenous absorbers and injected probes

• Depth & spatiotemporal resolution:

Acquire whole-body images of small animals, with in-plane resolution of 150 μm , in vivo and in real time




MSOT applications

Broad range of applications including cancer, cardiovascular, neuro and pharmacokinetic imaging

Cancer: Understanding tumors

- Measure probe accumulation, analyze distribution
- Determine oxygenation status
- Visualize and quantify (micro-)vasculature



Vascular perfusion
in tumor region

Cardio: Monitoring blood flow


- Visualize blood vessels in all body regions
- Quantify blood oxygenation in vessels and tissue
- Analyze atherosclerotic plaque deposition



Oxygenation
in heart region

Neuro: Tapping into the brain


- High resolution imaging of deep brain structures through intact skin and skull
- Determine probe distribution in the brain
- Characterize blood brain barrier integrity



Probe distribution
within the brain

Kinetics: Measuring the change

- Follow probe uptake and elimination in real time
- Analyze hepatic and renal function
- Kinetic profiling of multiple probes simultaneously



Probe accumulation and
clearance in kidneys

MSOT: SYSTEM OVERVIEW, IMAGING WORKFLOW, ACCESSORIES/OPTIONS

System overview



- Ergonomic and fully integrated design (no external components)
- Small footprint (< 1 sqm)
- System fully upgradeable



Workflow-oriented software design:

1. Acquisition
2. Reconstruction
3. Spectral processing
4. Visualization and analysis

Imaging workflow

Animal preparation



- Safe animal handling
- Repeatable fixed position (key for longitudinal studies)
- Easy mounting of animal in holder
- Integrated anesthesia supply

- Holder snapped into stage-driven sled
- Optimal animal positioning in three dimensions
- No limitations on regions of interest
- Imaging of entire animal in one session

Insertion into imaging chamber



Image acquisition



- x-sectional scan along the region of interest
- Fully automated image acquisition
- Access for catheter or thermometer
- Maintains animal body temperature

Accessories/options



Anesthesia system
Copyright © UNO BV 2012



Animal holder for imaging rats
(up to ~150 g)

Additional options include:

- Different laser wavelength ranges
- Advanced software functionality
- CCD camera for live animal monitoring
- Reference mouse anatomy atlas

MSOT: SCANNER FAMILY



Technical specifications	MSOT inSight 64		MSOT inVision 128	MSOT inVision 256-TF
Image acquisition				
Image rate (acquisition and live display)	10 fps (frames per second)			
Acquisition time, cross-sectional image	< 100 ms (single wavelength) / < 1s (multispectral)			
Spatial resolution	150 µm			
Lower limit of quantification (LLOQ), e.g. for:				
ICG	50 nM			
Gold nanoparticles	15 fM			
Key system components				
SpectraPULSE™ illumination system				
Wavelength spectrum (standard)	680-980 nm			
Pulse repetition rate	10 Hz			
Maximum pulse energy	80 mJ	100 mJ		
Wavelength tuning				
Tuning speed	< 100 ms			
Minimum step size	1 nm			
TomoARC™ detection system				
Curvature of detector array	172°	270°		
Center frequency / bandwidth	5 MHz / > 55 %			
Number of detector elements	64	128	256	
RapidSCAN™ data acquisition electronics				
Channels for simultaneous acquisition	64 (max: 512)	128 (max: 512)	256 (max: 512)	
Sampling rate	40 MSps (mega-samples per second)			
General technical specifications				
Scanner console - hardware	Intel Core i7, 5 TByte HDD, RAID, 24 GByte RAM, 1 GBit Ethernet, 24" TFT screen			
Scanner console - operating system	64-bit Windows 7			
Dimensions (width x depth x height)	89 x 95 x 165 cm			
Safety certification	CE marking (system compliant with IEC 61010-1, 60825-1, 61326-1)			
Laser classification	Class 1			

Software functionality


ViewMSOT™	Data management, image acquisition, image reconstruction, spectral unmixing, signal quantification, data/image export
FastPK™	Ultrafast multispectral data acquisition, PK analysis toolset, video creation and export
QuantumPRO™	Advanced toolset for signal quantification, component analysis, extraction of measurement reports
PulseCTRL™	Adjustment of pulse energy and repetition rate, advanced laser performance monitoring
CryoMOUSE™	Annotated cryoslice images for anatomical reference, male and female anatomy, synchronization with imaging session

MSOT: SYSTEM CONFIGURATION/ OPTIONS

System configuration/ options	MSOT inSight 64	MSOT inVision 128	MSOT inVision 256-TF
Motorized in-plane object translation (x+y axes)	<input type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>
Extended axial FOV (180 mm; including holder)	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Live animal monitoring (CCD camera inside imaging chamber)	<input type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>
Anesthesia flow monitoring	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Extended laser wavelength range (410-980 nm, 680-2.300 nm, 410-2.300 nm)	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Laser-integrated pulse energy meter	<input type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>
Laser-integrated spectrometer	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
FastPK™ analysis of fast kinetic processes¹	<input type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>
QuantumPRO™ advanced data analysis¹	<input type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>
PulseCTRL™ advanced laser control¹	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
CryoMOUSE™ reference anatomy²	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
Legend	<input checked="" type="radio"/>	Standard configuration	
	<input type="radio"/>	Optional configuration	

¹Requires laser-integrated pulse energy meter

²Requires live animal monitoring (CCD camera inside imaging chamber)



MSOT: PUBLICATIONS

Selection of academic publications summarizing the MSOT performance:

- X. Luís Deán-Ben and Daniel Razansky,
Adding fifth dimension to optoacoustic imaging: volumetric time-resolved spectrally enriched tomography,
Light: Science & Applications (2014) 3, e137; doi:10.1038/lisa.2014.18.
- Buehler A et al.,
Real-time handheld multispectral optoacoustic imaging,
Opt Lett. 2013 May 1;38(9):1404-6. doi: 10.1364/OL.38.001404.
- Stritzker J et al.,
Vaccinia Virus-mediated Melanin Production Allows MR and Optoacoustic Deep Tissue Imaging and Laser-induced Thermo-therapy of Cancer,
PNAS February 26, 2013 vol. 110 no. 9 3316-3320.
- Taruttis A et al.,
Multispectral Optoacoustic Tomography of Myocardial Infarction,
Photoacoustics (Vol. 1, Issue 1, March 2013). DOI: 10.1016/j.pacs.2012.11.001.
- Burton NC et al.,
Multispectral Opto-acoustic Tomography (MSOT) Brain Imaging and Characterization of Glioblastoma,
Neuroimage, 2012 Sep 28; pii: S1053-8119(12)00963-9.
- Herzog E et al.,
Optical Imaging of Cancer Heterogeneity with MSOT,
Radiology. 2012 May;263(2):461-8.
- Taruttis A et al.,
Fast Multispectral Optoacoustic Tomography (MSOT) for Dynamic Imaging of Pharmacokinetics and Biodistribution in Multiple Organs,
PLoS ONE 2012, 7(1):e30491.
- Razansky D et al.,
Volumetric Real-time Multispectral Optoacoustic Tomography (MSOT) of Biomarkers,
Nature Protocols 6, 1121-1129 (2011).
- Ntziachristos V and Razansky D,
Molecular Imaging by Means of Multispectral Optoacoustic Tomography (MSOT),
Chemical Reviews, 110(5), 2783-2794 (2010).