

### Anesthesia setup guide







3/16 id tubing



Oxygen hose





Table top anesthesia machine With flow meter



Vaporizer





F/air canister for evacuation



Rodent circuit set, shown With nose cone and diaphragm





manifold



Anesthetizing chamber



Regulator for O2 tank



#### Step one



+





Attach the oxygen regulator to the green oxygen hose





#### Step two



Attach the other end of the green oxygen hose to the bottom of the flow meter on the anesthesia machine





#### Step three







Mount the vaporizer on the anesthesia machine Using the screws provided, next attach tubing from The top of the flow meter to the inlet end cap of the Vaporizer as shown.







#### Step four



Connect the rodent circuit set to the outlet end cap of The vaporizer using a length of 3/16 id tubing as shown







### Final step








Lastly, connect the F/air canister to The rodent circuit set as shown. NOTE: when using the F/AIR canister it must be laying on its side





### Utilizing a manifold in your anesthesia setup

If you are planning on using a manifold to run multiple lines of anesthetic agent to various devices, simply set up as shown and change the following: 1) connect one side of the manifold to the outlet of the vaporizer instead of the rodent circuit set using either ¼" or 3/16" ID tubing, 2) connect the other side of the manifold to an F/AIR canister using the Tubing adapter, 3) connect anesthetizing box, rodent circuit set, gas anesthesia platform or other devices directly to the manifold with the 3/16 id tubing, 4) make sure the outlet on each device is connected to the side of the manifold which outlets to the F/AIR canister



If you have any questions about setting up your anesthesia system please contact:

Harvard Apparatus Technical Support at 1-800-272-2775