

# Fungi/Yeast Genomic DNA Isolation Kit


CAT. 27300, 27350

FOR THE RAPID PURIFICATION OF DNA FROM YEAST CELLS AND FUNGAL SPORES OR MYCELIUM




- ✓ Rapid spin column purification of genomic DNA from viable yeast cells, fungal spores or mycelium, and even Gram-positive bacteria
- ✓ Bead tubes (provided) allow for effective mechanical homogenization
- ✓ Purified DNA is of high quality and integrity and compatible with sensitive downstream applications such as PCR, qPCR, RFLP and more
- ✓ 96-Well format available for high throughput isolation

# Advantages You Will Bring to Your Lab



## EFFICIENT PURIFICATION OF DNA FROM DIFFERENT FUNGAL SPECIES


**Figure 1. DNA Isolation from Different Fungi Species and Yeast.** To demonstrate the purification of DNA from different fungal species, 30 mg of fungi were collected from plate cultures of *Pichia* sp., *Aspergillus niger*, *Cladosporium cladosporioides*, *Botrytis cinerea* and *Mucor racemosus*, and the DNA was extracted using Norgen's Fungi/Yeast Genomic DNA Isolation Kit. The bead system efficiently lysed the fungal cell walls with the provided Lysis Solution, and total DNA was eluted in 100  $\mu$ L. For analysis, 10  $\mu$ L from each elution was loaded in 1% 1xTAE agarose gel. Lane 1: Yeast (*Pichia* sp.), Lane 2: *Aspergillus Niger*; Lane 3: *Cladosporium cladosporioides*; Lane 4: *Botrytis cinerea*; Lane 5: *Mucor racemosus*; Lane M: Norgen's HighRanger 1kb DNA Ladder. The optional RNase treatment was not performed during the process.




## SUCCESSFULLY ISOLATE DNA FROM BOTH THE SPORES AND MYCELIUM.

**Figure 2. DNA can be Isolated from Spores and Mycelium.** Genomic DNA was extracted from the spores and mycelium of *Botrytis cinerea* using Norgen's Fungi/Yeast Genomic DNA Isolation Kit. For the spores,  $10^5$  spores/mL were collected from a plate culture by adding 3 mL of 0.9% NaCl to the fungi culture plate and gently shaking to extract the spores. For mycelium, 20 mg of mycelium was collected with sterilized forceps. DNA was isolated using Norgen's Fungi/Yeast Genomic DNA Isolation Kit with the provided Bead Tubes. For analysis, 10  $\mu$ L from each 100  $\mu$ L elution was loaded on a 1% 1xTAE agarose gel. Lane M: Norgen's HighRanger 1kb DNA Ladder. Optional RNase treatment was not performed during the extraction. As it can be seen, DNA could be successfully isolated from both the spores and mycelium.

## HIGH THROUGHPUT DNA ISOLATION WITHOUT COMPROMISING THE PURITY AND QUALITY



**Figure 3. High Yields of DNA Isolated from *Saccharomyces cerevisiae* overnight culture.** DNA was isolated from *S. cerevisiae* overnight culture using Norgen's Fungi/Yeast Genomic DNA Isolation 96-Well Kit. Following isolation of 10 samples, 10  $\mu$ L from each 100  $\mu$ L elution was loaded on 1% TAE agarose gel. Norgen's Fungi/Yeast Genomic DNA Isolation 96-Well Kit demonstrated a good and consistent DNA yield and integrity. Lane M: Norgen's HighRanger 1kb DNA Ladder.



**Figure 4. High quality of yeast DNA measured by NanoDrop.** DNA isolated from *Saccharomyces cerevisiae* overnight culture using Norgen's Fungi/Yeast Genomic DNA Isolation 96-Well Kit showed a high quality of DNA based on A260/280 and A260/230. This indicates that the isolated DNA from all samples is of a high quality and can be used in sensitive downstream applications.

## Kit Specifications

Description	Specifications
Column Binding Capacity	50 $\mu$ g
Maximum Column Loading Volume	650 $\mu$ L
Maximum Amount of Starting Material	
Fungi (wet weight)	50 mg
Yeast or Bacteria Culture	0.5mL - 1mL
Average Yield*	
Pichia sp. (yeast)	25 $\mu$ g
Botrytis cinerea	32 $\mu$ g
Fusarium sp.	42 $\mu$ g
Aspergillus niger	26 $\mu$ g
Mucor racemosus	15 $\mu$ g
Cladosporium cladosporioides	12 $\mu$ g
Penicillium sp.	40 $\mu$ g
Rhizopus oryzae	7 $\mu$ g
Alternaria tenuissima	30 $\mu$ g
Time to Complete 10 purifications	30 minutes

\* Yield will vary depending on the type of sample processed

## Select Publications

Bitterman, Roni, et al. **“Exhaled Breath Condensate Surveillance for Aspergillus in Acute Leukemia—a Pilot Trial.”** Open Forum Infectious Diseases, vol. 11, no. 10, Oct. 2024, p. ofae537. Silverchair, <https://doi.org/10.1093/ofid/ofae537>.

Duarte-Silva, Luciana, et al. **“Phenotypic and Molecular Characterization of Prototheca Wickerhamii from a Brazilian Case of Human Systemic Protothecosis.”** PLOS Neglected Tropical Diseases, vol. 18, no. 11, Nov. 2024, p. e0012602. PLoS Journals, <https://doi.org/10.1371/journal.pntd.0012602>.

García-Salazar, Eduardo, et al. **“Utility of Cand PCR in the Diagnosis of Vulvovaginal Candidiasis in Pregnant Women.”** Journal of Fungi, vol. 11, no. 1, 1, Jan. 2025, p. 5. www.mdpi.com, <https://doi.org/10.3390/jof11010005>.

Tran, Tra-My Thanh, et al. **“Isolation, Screening, and Identification of Thermophilic Lignocellulolytic Fungi from Agricultural Soil in Vietnam.”** IOP Conference Series Earth and Environmental Science, vol. 1399, no. 1, 1 Sept. 2024, pp. 012003–012003, <https://doi.org/10.1088/1755-1315/1399/1/012003>.

Öncel, Sibel, and Hilal Özkılınc. **“Discovering the Dynamics of Peach Fruit Mycobiome throughout Fruit Development Season by High-Throughput Sequencing.”** Scientific Reports, vol. 15, no. 1, Mar. 2025, p. 8969. www.nature.com, <https://doi.org/10.1038/s41598-025-93090-6>.

Preston, Michael et al. **“Expanded geographical distribution of Sarcosoma globosum (witches cauldron) in Canada.”** Acta Mycologica, vol. 59, 2024, pp. 1-12. doi:10.5586/am/195530.

Zenelt, Weronika, and Krzysztof Krawczyk. **“Insect-Derived Bacteria as Biocontrol Tool and a Potent Suppressor of Plant Pathogenic Fungi in Tomato Cultivation.”** Microbial Pathogenesis, vol. 198, Jan. 2025, p. 107158. ScienceDirect, <https://doi.org/10.1016/j.mic-path.2024.107158>.

## Ordering Information

Description	Preps	Cat. #
Fungi/Yeast Genomic DNA Isolation Kit	50 Prep	27300
Fungi/Yeast Genomic DNA Isolation Kit	2 x 96-Well plate	27350

## Related Products

Description	Prep Size	Cat. #
Bacterial Genomic DNA Isolation Kit	50 Prep	17900
Plant RNA/DNA Purification Kit	50 Prep	24400
Norgen Next Generation Sequencing Services	Call 1-866-Norgenb or Visit <a href="https://norgenbiotech.com">norgenbiotech.com</a>	



### Order Today!

Order Norgen's Fungi/Yeast Genomic DNA Isolation Kit today to help improve your research performance and workflow.

#### Three ways to order

Visit [norgenbiotech.com](https://norgenbiotech.com)  
Call at 1-866-NORGENB or 905-227-8848  
Email at [orders@norgenbiotech.com](mailto:orders@norgenbiotech.com)



For more information  
**SCAN HERE**

