

# Dry Bath

## HB120-S

### Dry bath

#### Features:

- A wide range of temperature control up to 120°C
- Timer function support is available
- Various types of optional blocks are available
- Overheating protection
- Lid for protection and heat preservation

#### Specifications

Specifications	HB120-S
Functions	Heating
Temperature range	Room temp.~120°C
Temperature control accuracy	± 0.5°C
Temperature uniformity	± 0.5°C
Max heating rate	5.5°C/min
Timer	1min~99h59min
Screen	LED
Overheating protection	140°C
Adapter block material	aluminum
Voltage, Frequency	100~120V/220~240V,50Hz/60Hz
Power	160W
Dimension[WxDxH ]	95x150x50mm
Weight	3kg

#### Cat.No.

#### Block types

#### Capacity

#### Diameter Depth

18900218



0.2mLx54

6.4mm  
15mm

18900219



0.5mLx40

8.2mm  
26.4mm

18900253

(1.5ml)

18900220

(2ml)



1.5/2mLx40

11mm  
34.7mm

18900221



5/15mLx28

17mm  
48mm

18900222



50mLx8

29.8mm  
48mm

18900224



0.2mL+0.5mL+1.5/2mL  
18+18+18

-

18900223



96 or 384 well plate

-



Adapter block selection



Overheating protection



Timer



LED display

