

Orbital and Linear Shaker

SK-O180-E/SK-L180-E





Feature:

- Compact, orbital shakers with ideal shaking motion, for a maximum shaking weight of 3kg
- Compact space-saving design fits easily in the incubator and refrigerator
- Continuous operation, speed range of 0-200 rpm
- Stable operation

Model	SK-O180-E/SK-L180-E
Voltage[VAC]	100-240
Frequency[HZ]	50/60
Shaking motion	Orbital/Linear
Orbital diameter[mm]	20
Max. load capacity (with platform)[Kg]	3
Motor rating input[W]	16
Motor rating output[W]	10
Speed range[rpm]	0-200
Dimensions[W×H×D mm]	330X270X110
Weight[kg]	3.1
Permissible ambient temperature[°C]	5-40
Permissible relative humidity	80%
Protection class acc. DIN EN60529	IP21

Ordering Information

P/N	Descriptions
831112010000	SK-O180-E, US plug, 100V-240V, 50Hz/60Hz
831112110000	SK-O180-E, CN plug, 100V-240V, 50Hz/60Hz
831112210000	SK-O180-E, Euro plug, 100V-240V, 50Hz/60Hz
831112310000	SK-O180-E, UK plug, 100V-240V, 50Hz/60Hz
832112010000	SK-L180-E, US plug, 100V-240V, 50Hz/60Hz
832112110000	SK-L180-E, CN plug, 100V-240V, 50Hz/60Hz
832112210000	SK-L180-E, Euro plug, 100V-240V, 50Hz/60Hz
832112310000	SK-L180-E, UK plug, 100V-240V, 50Hz/60Hz