



90 x 14 mm Petri Dish

Made in polystyrene. Supplied in groups of 20 units, packaged in heat sealed bags.
Code **200200** is **aseptic**. Code **200209** is **sterile by radiation**.
Suitable for automatic filling machines.

code	description	sterile	case quantity	case weight	case volume	cases per pallet
------	-------------	---------	---------------	-------------	-------------	------------------

With three vents

200200	Ø 90 x 14 mm	no	25 x 20	7.32	0.071	28
200209	Ø 90 x 14 mm	STERILE R	25 x 20	7.39	0.071	28

Not vented dishes, for the cultivation of anaerobic

200200.4	Ø 90 x 14 mm	no	25 x 20	6.90	0.071	28
-----------------	--------------	----	---------	------	-------	----



90 x 14 mm Petri Dish, two compartments

Made in polystyrene. **Aseptic**.
With three vents.
Supplied in groups of 20 units, packaged in heat sealed bags.
Suitable for automatic filling machines.

code	description	case quantity	case weight	case volume	cases per pallet
200202	Ø 90 x 14 mm 2 compartments	25 x 20	7.75	0.072	28



90 x 14 mm Petri Dish, three compartments

Made in polystyrene. **Aseptic**.
With three vents.
Supplied in groups of 20 units, packaged in heat sealed bags.
Suitable for automatic filling machines.

code	description	case quantity	case weight	case volume	cases per pallet
200203	Ø 90 x 14 mm 3 compartments	25 x 20	7.82	0.070	28



90 x 15 mm Petri Dish, four compartments

Made in polystyrene. **Aseptic**.
Vented.
Supplied in groups of 24 units, packaged in heat sealed bags.
Suitable for automatic filling machines.

code	description	case quantity	case weight	case volume	cases per pallet
200210	Ø 90 x 15 mm 4 compartments	24 x 25	9.20	0.077	28



90 x 25 mm Petri Dish, four compartments

Made in polystyrene. **Aseptic**.
Vented.
Supplied in groups of 24 units, packaged in heat sealed bags.
Suitable for automatic filling machines.

code	description	case quantity	case weight	case volume	cases per pallet
200215	Ø 90 x 25 mm 4 compartments	24 x 13	8.63	0.076	24

Petri Dish 90 x 15 mm with internal cross

Made in polystyrene. **Aseptic.**

With four vents. Stackable.

With a mark on a side that helps orientation and a writing area on the bottom.

Supplied in groups of 25 units, packaged in heat sealed bags.

Suitable for automatic filling machines.

Internal cross guarantees the totally flat base.

code	description	case quantity	case weight	case volume	cases per pallet
200200.5	Ø 90 x 15 mm internal cross	24 x 25	7.50	0.076	20



140 x 20 mm Petri Dish

Made in polystyrene. With three vents.

Supplied in groups of 10 units, packaged in heat sealed bags.

Suitable for automatic filling machines.

Code **200214** is **aseptic**.

Code **200219** is **sterile by radiation**.

code	description	sterile	case quantity	case weight	case volume	cases per pallet
200214	Ø 140 x 20 mm	no	11 x 15	7.40	0.074	30
200219	Ø 140 x 20 mm	STERILE R	11 x 15	7.25	0.072	30



120 x 120 mm squared Petri Dish

Made in polystyrene. **Aseptic.**

With four vents.

Supplied in groups of 10 units, packaged in heat sealed bags.

Suitable for automatic filling machines.

code	description	case quantity	case weight	case volume	cases per pallet
200204	120 x 120 mm	24 x 10	11.18	0.072	28



55 x 14 mm Petri Dish

Made in polystyrene. With four vents.

Supplied in groups of 15 units, packaged in heat sealed bags.

Suitable for automatic filling machines.

Code **200201** is **aseptic**.

Code **200201.B** is **sterile by radiation**.

code	description	sterile	case quantity	case weight	case volume	cases per pallet
200201	Ø 55 x 14 mm	no	80 x 15	8.40	0.068	30
200201.B	Ø 55 x 14 mm	STERILE R	80 x 15	8.40	0.068	30

