

Blood Bank Refrigerators

HXC-358B



Ergonomic Design



Famous Compressor



Stainless Steel Drawer



Condensation Free Design



Intelligent Alarm



Infinite Divided Stainless Steel Drawer

Stainless steel drawer can be unlimited separated based on the demand



Three Types of Alarm

Audible buzzer, visible flashing light, remote alarm



Five Alarm Functions

High/low temperature, power failure, sensor error, low battery, door ajar

Specifications

HXC-358B

Technical Data

| | |
|---------------|--------------------|
| Type | Drawer-Type |
| Climate Class | ST |
| Cooling Type | Forced air cooling |
| Defrost Mode | Auto |
| Refrigerant | CFC-Free |
| Noise((dB(A)) | 43 |

Performance

| | |
|-------------------------|-------|
| Temperature Range(°C) | 2~6 |
| Ambient Temperature(°C) | 10~38 |

Control

| | |
|------------|----------------|
| Controller | Microprocessor |
| Display | LED |

Electrical Data

| | |
|-----------------------|---------------|
| Power Supply(V/Hz) | 220-240/50/60 |
| Power(W) | 460W |
| Electrical Current(A) | 3 |

Dimensions

| | |
|--|-------------------------------------|
| Capacity(L/Cu.Ft) | 358/12.7 |
| Blood Storage Capacity (450ml blood bags) | 200 |
| Net/Gross Weight(approx) | 165/181(kg) 363.8/399.0(lbs) |
| Interior Dimension(W*D*H) | 620*490*1160(mm) 24.4*19.3*45.7(in) |
| Exterior Dimension(W*D*H) | 720*690*1730(mm) 28.3*27.2*68.1(in) |
| Packing Dimension(W*D*H) | 795*805*1870(mm) 31.3*31.7*73.6(in) |
| Container Load(20'/40'/40'H) | 14/30/30 |



Functions

| | |
|----------------------|---|
| High/Low Temperature | Y |
| Power Failure | Y |
| Sensor Error | Y |
| Low Battery | Y |
| Door Ajar | Y |
| Remote Alarm | Y |

Accessories

| | |
|---------------------|----------|
| Caster | Y |
| Foot | Y |
| Test Hole | Y |
| Shelves /Basket | - |
| Drawers/Inner doors | 5/- |
| USB Interface | Optional |
| Temp Recorder | Y |

Others

| | |
|-------------|----|
| Certificate | CE |
|-------------|----|

Haier Biomedical UK

Address: 2 Birch House, Brambleside, Bellbrook
Industrial Estate, Uckfield, East Sussex,
TN22 1QQ, United Kingdom
TEL: +44(0) 1825 744 696 (Head Office)
www.haierbiomedical.co.uk



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.



Scan 2D barcode with your mobile device to learn more information about Haier Biomedical