

- Conforms to the requirement of EN 285 European Standard for Large Steam Sterilizers
- Chamber and STU capacity:
NC 300: 300 liters, 4 pcs. STU
NC 430/430D: 430 liters, 6 pcs. STU
NC 570/570D: 570 liters, 8 pcs. STU
NC 710/710D: 710 liters, 10 pcs. STU
- Ideal design for the steam sterilization of packed or unpacked surgical and dental instruments, operation linen, glass, plastic, rubber and silicon materials, infusion liquids, microbiological cultures and medical waste
- Chamber, steam jacket, door and steam generator made of 316L stainless steel and conform to 2014/68/EU Pressure Equipment Directive
- Pneumatically controlled door opening with safety system
- Steam driven one piece chamber gasket made of silicon
- High capacity, low noise water ring vacuum pump for efficient vacuum
- High quality valves, sensors and transmitters for a problem free and safe operation
- Pneumatic valves in the steam lines
- Innovative AIR – COOLING system for low water and energy consumption
- N-Wise™ microprocessor control system with 800x480 pixels 7" colourful touch screen
- Wide screen allows monitoring all related information regarding sterilization cycle
- 7 preset programs, 50 custom made programs
- Two test programs: Bowie & Dick and vacuum leakage
- Memory for last 1.000 cycles
- Independent control system besides the main control system to record the data related to every cycle
- Comprehensive self-diagnostic system for possible system malfunctions
- Service menu to detect the problem in case of a failure
- Thermal printer offered as standard
- Standard NuveCom™ communication unit
- Sending e-mails up to five e-mail addresses with the details of failure
- USB port to store data on memory stick
- World-wide access via Ethernet port with NuveCloser™
- Equipped with Wi-Fi connection for N-Mobile™
- Optional AlerText™ GSM alarm module
- Password protected menu access
- Standard RS 485 connection for PC
- Double door versions with 4,7" display on the clean side:
NC 430D, NC 570D and NC 710D

**NC 300/430/570/710
NC 430D/570D/710D
STEAM STERILIZERS
POWERED BY N-Wise™**

