

96-Well Magnetic Plate

Product # 75600

Product Insert

Product Description

Norgen's 96-Well Magnetic Plate is designed for easy manual nucleic acid isolations, and other magnetic-bead based purifications including cleanup and library preparation. The plate has 4 corner tops with strong permanent magnets to secure 96-well plates for rapid, fast separations. The 96-Well Magnetic Plate can sit flat on a lab bench for manual isolations.

Using the 96-Well Magnetic Plate

The 96-Well Magnetic Plate is intended to be used with magnetic bead-based purification kits. It is recommended to use Norgen's magnetic bead-based purification kits. Norgen's 96-Well Magnetic Plate is compatible with most of the 96-Well plates.

Kit Components

Component	Quantity
96-Well Magnetic Plate	1
Product Insert	1

Advantages

- Strong and permanent magnets with high magnetic field
- Fast, rapid separations
- Can hold 96 Well Deep well plate.
- Can accommodate plates with volumes from 50 μ L to 2 mL
- Ergonomically designed
- Easy to use
- Lightweight and durable

Storage Conditions and Product Stability

Keep away from strong and external magnets. Store at room temperature in dry conditions free from dust and other contaminants.

Precautions and Disclaimers

This kit is designed for research purposes only. It is not intended for human or diagnostic use. Ensure that a suitable lab coat, disposable gloves and protective goggles are worn when working with chemicals. Strong magnets in the plates can interfere with medical devices such as pacemakers. Individuals with these implants should stay away and maintain a distance from the plates.

Related Products	Product #
Plant DNA Isolation 96-Well Kit (Magnetic Bead System)	62400
Cells and Tissue DNA Isolation 96-Well Kit (Magnetic Bead System)	62500
Blood DNA Isolation 96-Well Kit (Magnetic Bead System)	62600
PCR and Sequencing Reaction Clean-Up 96-Well Kit (Magnetic Bead System)	62700
Soil DNA Isolation 96-Well Kit (Magnetic Bead System)	62800
Plasmid MiniPrep 96-Well Kit (Magnetic Bead System)	63000
Stool DNA Isolation 96-Well Kit (Magnetic Bead System)	63100
Saliva DNA Isolation 96-Well Kit (Magnetic Bead System)	RU62900
Adhesive Tape	28394

Technical Support

Contact our Technical Support Team between the hours of 9:00 and 5:30 (Eastern Standard Time) at (905) 227-8848 or Toll Free at 1-866-667-4362.

Technical support can also be obtained from our website (www.norgenbiotek.com) or through email at support@norgenbiotek.com.

3430 Schmon Parkway, Thorold, ON Canada L2V 4Y6
Phone: (905) 227-8848
Fax: (905) 227-1061
Toll Free in North America: 1-866-667-4362